

NEW

Introducing

ZEON[®]

ENDEAVOUR[™]

LIGHTWEIGHT COMFORT AND POWER COMBINED



Actual Size

Orascoptic's lightest & smallest light system ever!

For over 30 years, Orascoptic has been a leader in lighting innovation. Applying advanced technology, Orascoptic has created the Endeavour[™], a high resolution headlight system that incorporates today's best technology at an affordable price.



Compact and lightweight for comfort

Designed with comfort in mind, the compact light and battery unit allow you to focus on your patients.



Innovative "Gumdrop" technology

The miniature light delivers a high quality focused beam for better visibility during treatment and diagnosis.



Advanced Light-Touch controls

Allows you to easily and quickly adjust light settings without reaching into a coat or pants pocket.



To see a demonstration of Endeavour[™]
or to buy online go to:
orascoptic.com/endeavour

Orascoptic

Sybron Dental Specialties

Superior Visualization[™]

3225 Deming Way, Suite 190 • Middleton, WI 53562

Toll Free: 800.369.3698 • 262.831.2555

Orascoptic's lightest & smallest light system ever!

ZEON[®] ENDEAVOUR[™]

LIGHTWEIGHT COMFORT AND POWER COMBINED



*Four-year warranty on the LED and cable. One-year warranty on the battery unit.



Battery power level display indicates how long before recharge is necessary, **ensuring procedures can be completed.**

Capacitive Touch controls for easy on/off and light adjustment - **no fumbling, buttons or switches.**

Small, lightweight battery pack for **maximum comfort.** Place in a pocket or on an armband.

High resolution illumination creates **better visibility for improved procedural effectiveness.**

Low heat output means **safety for both patient and practitioner.**

The smallest headlight to **deliver at least 4,000 foot candles of white light.**

The **high quality** focused spot provides **superior resolution** for more efficient treatment.

Product Specifications

HEADLIGHT SPECIFICATIONS

Diameter:	0.59 in (15 mm)
Length:	0.63 in (16 mm)
Weight:	0.26 oz (7.37 grams)
Light Intensity:	2000/3000/4000 ft cd of light output
Intensity Control:	Capacitive Touch up to 4000 ft cd
Focus:	2.25" (5.7 cm) spotlight at 12" (30.5 cm) focused beam
Composite Filter:	Curing filter included

BATTERY SPECIFICATIONS

Dimensions:	3.5" (8.9 cm) x 2" (5.1 cm) x 0.8" (2 cm)
Weight:	4.6 oz (115 grams)
Est. Battery Life:	5 to 10 hours depending on light setting
Est. Charging Time:	3 hours (can be used while recharging)
Type:	Lithion Ion (the highest capacity available)

Orascoptic
Sybron Dental Specialties



Superior Visualization[™]

3225 Deming Way, Suite 190 • Middleton, WI 53562

Toll Free: 800.369.3698 • 262.831.2555

orascoptic.com/endeavour