

INTRODUCING

GEM

QUALITY. LIGHT-WEIGHT. AFFORDABLE.



All Orascope Loupes are:

- Customized to your facial geometry, working distance, and angle of declination
- Hypo allergenic
- Corrosion resistant

Loupes Package Includes:

- Side shields
- Head strap
- Sleek black storage case
- Prescription lenses (if applicable)
- Protective lens caps

Designer Victory Frame

Material: Titanium
Sizes: Small, Medium, Large
Unlimited Prescriptions



Slate

Wine

Magnification: Class II
Field Width: 5" (13 cm)
Field Depth: 6" (15 cm)
Weight: 1.9 oz. (54 grams)
Working Distance:
 Short: 10" – 16" (26 – 41 cm)
 Regular: 14" – 20" (36 – 51 cm)

Orascope's new **GEM** is the lightest flip-up loupe on the market. Finally, dental professionals who prefer the convenience of flip-up magnification can practice all day in total comfort.

Renowned for **quality craftsmanship**, Orascope loupes feature a wide and deep viewing field and a limited lifetime warranty.

In addition to being **exceptionally lightweight**, this little **GEM** incorporates new optical materials that make it a **highly affordable** addition to the Orascope portfolio of visualization products.



Superior Visualization™

Kerr Corporation

3225 Deming Way, Suite 190, Middleton, WI 53562 800.369.3698 • T: 608-831-2555 www.orascope.com

©2010 Kerr Corporation

960836 Rev. B 4/10

