

# Orasoptic DK™

SEE WHAT YOU'VE BEEN MISSING



**Orasoptic**

*An innovative diagnostic kit for dental professionals.*

# Orasoptic DK™

PROFESSIONAL DENTAL DIAGNOSTIC KITS



## Endo Fiber Lite

use with Fiber Lite Instrument

This disposable fiber aids in detection of caries with high definition. Root fractures are also easily detected by placing the fiber in the gingival sulcus.



## Transilluminator - 3mm

High quality glass light guide which allows doctors to better visualize caries on anterior and posterior teeth.

## Barrier Sleeves *not shown*

Disposable clear sleeves that fit over the transilluminator to help prevent cross-contamination.

## Autoclavable

All instruments are autoclavable. The LED Light Source and Endo/Perio Fibers themselves are not.

Orasoptic DK™ diagnostic instruments allow dental professionals to detect caries, root recession, and the oral cavity in general. The PRO kit includes the BASIC, ENDO, and PERIO kits have been designed for specific needs. Choose a kit to your specific needs.

## Orasoptic DK BASIC

Storage Case  
LED Light Source  
Transilluminator - 3mm  
Barrier Sleeves (250/bx)

## Orasoptic DK ENDO

Storage Case  
LED Light Source  
Fiber Lite Instrument  
(2) 25pc Packs of Endo Fibers



# See what you've been missing...

professionals to better visualize caries, root fractures, periodontal  
des the entire portfolio of DK instruments, while the  
specific dental applications. You can also custom-configure

## Orascope DK PERIO

Storage Case  
LED Light Source  
Fiber Lite Instrument  
(2) 25pc Packs of Perio Fibers

## Orascope DK PRO

Storage Case  
LED Light Source  
Transilluminator - 3mm  
Lighted Mirror  
Fiber Lite Instrument  
(1) 25pc Pack of Endo Fibers  
(1) 25pc Pack of Perio Fibers  
Barrier Sleeves (250/bx)



### Perio Fiber Lite

use with Fiber Lite Instrument

This disposable fiber helps measure periodontal recession. There are tick marks on the fiber at 3, 5, and 7 mm which, when lit, are easy to see.

### LED Light Source

Protected by an aluminum housing, the LED Light Source emits a beam of light that is used by all the Orascope DK™ Kit accessories.

### Lighted Mirror

Provides additional light for hard to reach areas.

### Fiber Lite Instrument

Used with Endo & Perio Fiber Lites.



**ORASCOPTIC DK™ PRE-ASSEMBLED KITS**

**ORASCOPTIC DK™ INDIVIDUAL ITEMS**

**Orascope DK BASIC**

**P/N 910832-1**

- Storage Case
- LED Light Source
- Transilluminator - 3mm
- Barrier Sleeves (250/bx)

**Orascope DK ENDO**

**P/N 910832-3**

- Storage Case
- LED Light Source
- Fiber Lite Instrument
- (2) 25pc Packs of Endo Fibers

**Orascope DK PERIO**

**P/N 910832-2**

- Storage Case
- LED Light Source
- Fiber Lite Instrument
- (2) 25pc Packs of Perio Fibers

**Orascope DK PRO**

**P/N 910832-4**

- Storage Case
- LED Light Source
- Transilluminator - 3mm
- Lighted Mirror
- Fiber Lite Instrument
- (1) 25pc Pack of Endo Fibers
- (1) 25pc Pack of Perio Fibers
- Barrier Sleeves (250/bx)

ITEM	P/N
LED Light Source	921755
DK™ Kit Case	921758
Lighted Mirror	921767
Transilluminator	921768
Fiber Lite Instrument	921769
Perio Fiber 25pk	921770
Endo Fiber 25pk	921771
Barrier Sleeve 250/bx	921772
Battery 3pk	921759



Call today to  
schedule a product  
demonstration:  
**800.369.3698**

**Kerr Corporation**  
3225 Deming Way, Suite 190  
Middleton, WI 53562  
T: 608.831.2555  
[www.orascope.com](http://www.orascope.com)



S u p e r i o r V i s u a l i z a t i o n ™