

BodyGuard[®] LX

ERGONOMIC SEATING SYSTEMS

Relieve physical stress due to:

- **Leaning forward** – rotates the pelvis backward reversing the normal forward bend in the lower spine.
 - This creates excessive pressure on the spinal discs and the body becomes unstable.
 - Uneven body weight distribution produces pressure under the thighs, impairing leg circulation.
- **Arms held away from the body** – places excessive stress on the musculature of the upper arms and shoulders.
- **Bending neck forward to gain access to the posteriors** – affects the musculature of the upper back and neck.

These are the critical factors addressed in the design of the **BodyGuard[®] LX** seating systems.

Extra cushioning and lumbar support.

Back tilts perfectly into the lower back for maximum support.

Luxurious. Stitched Ultra-Leather fabric. **Plush** Visco-Elastic cushioning.

Multiple arm adjustments.

Dual armrests move five directions to conform with the natural angle of clinician's arm, removing stress from back, neck and shoulders.



PROPER POSTURE means less stress to joints and muscles.

Fits all body sizes.

Three arm mounting positions and three height lifts fit all body sizes.

Seat tilt. Allows for continuous circulation in lower body during long operating procedures.



Operator's Chair & Armless Stool:

Short height	Regular height	Tall height
14" - 19" (36cm - 48cm)	17" - 23" (43cm - 58cm)	20.5" - 28.5" (52cm - 72cm)

Tall height Operator's Chair or Stool are available with or without a footring.

Assistant's Chair:

Seat height 20.5" - 28.5" (52cm - 72cm)
Left handed dentist: Right arm mount
Right handed dentist: Left arm mount



SINGLE COLOR
BodyGuard® LX chair



SIGNATURE DUAL-COLOR
BodyGuard® LX chair

The BodyGuard® LX line employs luxurious Ultra-Leather® fabric, soft Visco-Elastic® cushioning, and an upscale stitched pattern on both the seat and back rest.



(Over 40 custom colors are also available and subject to longer lead times)



Superior Visualization™