

INTRODUCING

# GEM

QUALITY. LIGHT-WEIGHT. AFFORDABLE.



**All Orascope Loupes are:**

- Customized to your facial geometry, working distance, and angle of declination
- Hypo allergenic
- Corrosion resistant

**Loupes Package Includes:**

- Side shields
- Head strap
- Sleek black storage case
- Prescription lenses (if applicable)
- Protective lens caps

**Designer Victory Frame**

**Material:** Titanium  
**Sizes:** Small, Medium, Large  
Unlimited Prescriptions



Slate

Wine

**Magnification:** Class II  
**Field Width:** 5" (13 cm)  
**Field Depth:** 6" (15 cm)  
**Weight:** 1.9 oz. (54 grams)  
**Working Distance:**  
 Short: 10" – 16" (26 – 41 cm)  
 Regular: 14" – 20" (36 – 51 cm)

Orascope's new **GEM** is the lightest flip-up loupe on the market. Finally, dental professionals who prefer the convenience of flip-up magnification can practice all day in total comfort.

Renowned for **quality craftsmanship**, Orascope loupes feature a wide and deep viewing field and a limited lifetime warranty.

In addition to being **exceptionally lightweight**, this little **GEM** incorporates new optical materials that make it a **highly affordable** addition to the Orascope portfolio of visualization products.



Superior Visualization™

Kerr Corporation 3225 Deming Way, Suite 190, Middleton, WI 53562 800.369.3698 www.orascope.com

©2009 Kerr Corporation

960836 Rev. A 10/09

